

# Shoprider

## 889D

ASK FOR THE BRAND SHOPRIDER



**SHOPRIDER**<sup>®</sup>  
*mobility*

[www.shoprider.com.au](http://www.shoprider.com.au)

RELIABLE MOBILITY GUARANTEED

# Shoprider

## 889D

**SHOPRIDER®**  
www.shoprider.com.au *mobility*

The SHOPRIDER™ 889D is scooter designed to assist operators who need a larger width scooter than is normally available. The scooter is well featured, it offers exceptional seating capabilities that includes a sliding seat for added flexibility. The 889 D has great features as standard including front and rear suspension, adjustable steering tiller, rear indicators, front light, sliding bench seat and fold up arm rests. The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

**THIS SCOOTER IS STRICTLY SINGLE PERSON USE ONLY ON AUSTRALIAN ROAD NETWORKS**

### SPECIFICATIONS

### 889D

Length	1480mm
Width	910mm
Height	980mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	12 Degrees
*Maximum Load <u>Single user</u>	225Kgs
Charger	8Amp
Motor Size	2.0Hp
Maximum Speed	12Kph
Battery capacity min	2 x 75 Amp
Maximum range#	up to 26km
Wheel Size	330mm
Turning Radius	2550 mm
Weight (with battery)	159kg


### FEATURES:

Wide bench Seat  
 Loop Steering Tiller  
 Adjustable Steering Tiller  
 Adjustable width arm rests  
 Sliding Seat  
 Front / Rear suspension  
 Battery Power indicator  
 Head Lights  
 Rear Lights & Indicators  
 Hazard Lights  
 Shopping Basket  
**Comprehensive Warranty**

**WARNING:** Do not exceed 10kph on public roads or paths. **SINGLE PERSON USE ONLY**  
**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

# Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



\*(Single user)

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

**SHOPRIDER®**  
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633  
 Or visit our Web Site at www.shoprider.com.au